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## English

### INSTRUCTIONS FOR USE AND PRESERVATION

Dear customer,

Thank you for your trust in this Orliman product. You have acquired a quality product with a high medicinal value. Please, read the instructions carefully. If you have any doubt, get in contact with your doctor, your orthopaedist specialist or our customer service department. Orliman appreciates your choice and wishes you a speedy recovery.

#### REGULATION

These products satisfy the European Medical Device Directive 93/42/CEE (RD. 1591/2009). A Risk Analysis has been carried out (UNE EN ISO 14971) minimising all the existing risks. The trials have been carried out in accordance with the European regulation UNE-EN ISO 22523 of Prosthesis and Orthesis.

#### SUITABLE FOR

Due to the variety of products available in the Thermo-med One Size range, the most important uses have been detailed. If you were to acquire a Thermo-med One Size not included in the following list, please consult your issuing establishment or your doctor's orders.

All the orthoses are suitable as a method of protection during sport activities, as well as for Thermo compression, stabilisation, protection, reduction of the oedemas, inflammation processes etc. Depending on the product:

**Ankle support:** Medio-lateral stabilisation, immobilisation, postural control.

**Knee support:** Centering the kneecap, prevention, protection for sport and work activities.

**Back support:** Lumbar pain, rheumatic pain and for activities that require continual strain.

**Elbow support:** Prevention of sport and work injuries, epicondylitis.

**Wrist support:** Stabilisation of the lower radio cubital joint, prevention, protection during sport and work activities.

#### FITTING INSTRUCTIONS

**THE PRODUCT SHOULD BE ADJUSTED BY YOUR ORTHOPAEDIC TECHNICIAN OR IN ACCORDANCE WITH DOCTOR'S ORDERS.**

In order to obtain a greater degree of therapeutic effectiveness in the different pathologies and to prolong the life and utility of the product, it is fundamental that the correct size adequate for each patient or user is chosen (enclosed in the packet is a diagram with all the sizes and the equivalent in centimetres). An excess compression can produce intolerance, therefore we recommend that the compression be regulated to an optimum degree.

#### To fit observe the following aspects:

Depending on the product:

**Ankle support:** Place the heel in the circular aperture. Tighten the narrowest pin of the foot from the sole towards the instep. Subsequently tighten the widest pin on the ankle.

**Knee support:** Place the centre of the opening over the kneecap, then tighten the thigh and calf pins on the posterior side of the leg (or knee).

**Back support:** Place the support on the back, 3cm from the gluteal fold, then proceed with the tightening and fastening.

**Elbow support:** Insert the arm through the pin to avoid any movement, centering over the elbow (perforated area). Tighten the arm and forearm pins.

**Wrist support (4607):** Insert the thumb through the interior ring, wrap the hand inserting it once again through the thumb tip, surround the wrist tightening the velcro fastener.

**Wrist support (4600):** Loosen the securement strap, allowing the maximum opening of the wrist guard, then introduce the hand and place it at the height of the wrist. Tighten the adjustment strap until the required securement is achieved. Finally fasten the strap using the velcro system.

**Cuff:** should be placed around the proximal forearm in pronation and the muscles relaxed about 3 cm below the epicondyle, allowing full flexion of the elbow to 90 degrees.

### PRECAUTIONS

The constructive material is inflammable.

Do not expose these products to situations where they could ignite. If this were the case remove them quickly and use the necessary measures to put it out. The material used is hypoallergenic, however we cannot guarantee 100% that in some determined cases it will not produce allergies, if this occurs, remove the product and see the doctor that prescribed it.

In case of minor irritations due to sweat, we recommend the use of a cotton sock, to separate the skin from contact with the material. For small inconvenience caused by sweat, we recommend using an interface to separate the skin from contact with the tissue. If rash, irritation or swelling to remove the product and consult a doctor or prosthetist. Contraindicated in open scars with swelling, redness and accumulation heat.

### RECOMMENDATION- WARNINGS

These products are to be used only by those which have been medically prescribed to do so. Persons unconnected to such prescriptions should not use them. When the product is not being used, store in its original packaging. In order to dispose of the packaging and the product, strictly comply with the legal norms of your community.

The use of the orthesis is conditioned to the recommendations of the prescribing doctor, therefore it must not be used for any other means than those prescribed.

For the validity of the guarantee, the issuing establishment should comply with this sheet of instructions.

### MANUFACTURE - CHARACTERISTICS


The material features used in production have been tested and approved and meet all the European quality requirements.

All the products are produced with high quality materials and offer an unbeatable comfort and quality of use. All products offer restraint, stability and compression for the optimum treatment of the pathologies for which they have been designed.

### RECOMMENDATIONS FOR PRESERVATION AND CLEANING

Stick the velcros to each other (if the orthesis has any), wash periodically by hand with a damp cloth and a neutral soap.

To dry the orthesis, use a dry towel to absorb most of the dampness and leave to dry at room temperature. Do not hang up or iron and do not expose to direct heat sources such as stoves, heaters, radiators, direct sun light etc. During use or during washing do not use alcohols, ointments or dissolvent liquids. If the orthesis is not properly dried any detergent residues could irritate the skin and deteriorate the product.

### GUARANTEES

ORLIMAN, S.L.U. Guarantees all its products, as long as these have not been manipulated or altered from their original state. Does not guarantee those products in which as a result of misuse, deficiencies or breakage of any kind, its original state has been altered. If you observe any deficiency or anomaly, communicate this immediately to the establishment from which it was obtained in order for it to be changed.

ORLIMAN, S.L.U. Appreciates your choice and wishes you a speedy recovery.

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