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English
INSTRUCTIONS FOR USE AND PRESERVATION

Dear customer,

Thank you for your trust in this Orliman product. You have acquired a quality product with a high medicinal value. Please, read the instructions carefully. If you have any doubt, get in contact with your doctor, your orthopaedist specialist or our customer service department. Orliman appreciates your choice and wishes you a speedy recovery.

REGULATION

These products satisfy the European Medical Device Directive 93/42/CEE (RD. 1591/2009). A Risk Analysis has been carried out (UNE EN ISO 14971) minimising all the existing risks. The trials have been carried out in accordance with the European regulation UNE-EN ISO 22523 of Prosthesis and Orthesis.

SUITABLE FOR

Instability of the knee, post-operative and post-trauma rehabilitation and relative immobilisation of the knee joint in situations that require a perfect control.

FITTING INSTRUCTIONS

THE PRODUCT SHOULD BE ADJUSTED BY YOUR ORTHOPAEDIC TECHNICIAN OR IN ACCORDANCE WITH DOCTOR'S ORDERS.

In order to obtain a greater degree of therapeutic effectiveness in the different pathologies and to prolong the life and utility of the product, it is fundamental that the correct size adequate for each patient or user is chosen (enclosed in the packet is a diagram with all the sizes and the equivalent in centimetres). An excess compression can produce intolerance, therefore we recommend that the compression be regulated to an optimum degree.

To fit observe the following aspects:

- 1-Open the fastening clips of the strap (a) and unfasten the buckles (b).
 - 2-Separate the hinge bars, extend the knee brace, and put on the device with the knee centred between the hinges. Orient the knee so that the hinges are pointing in the indicated direction, with the calf pads (the smaller ones) towards the feet.
 - 3-Put on the hinge bars laterally and medially with respect to the leg, and centre the hinge at the knee joint.
 - 4-Attach the two straps closest to the knee, without tightening.
 - 5-Loosen the friction clips. For correct adjustment, slide the upper and lower telescopic hinge bars in order to accommodate the leg along the entire length. Close the friction clips. The length indicators of the hinge bars allow verifying the concordance of the selected lengths on the thigh and on the calf.
 - 6-Attach the two remaining straps, without tightening.
 - 7-Pull on the straps to tighten them until there is no gap behind the leg. Be careful to maintain the lateral and medial positions of the hinge bars.
 - 8-Close the friction clips.
 - 9-Pull on the straps in order to tighten them through the buckles. Be careful to maintain the lateral and medial positions of the hinge bars.
 - 10-Fasten the ends of the straps. Use the Y-tongues at the ends of the straps in order to fasten them. The straps may have to be folded in order to shorten them before introducing them in the Y-tongues.
- Adjust the movement range control hinge:
- 11-The values of the extension limits can be selected between 10° of hyperextension and 30° by pulling the tongue out and sliding it to the desired position.
 - 12-The values of the flexion limits can be selected between 10° of hyperextension and 120°.
 - 13-The hinge can be fastened in one position by sliding the quick fastening button to the blocking position at any of the 5 positions: 10° of hyperextension, 0° (neutral), 10°, 20° and 30°.

For the model ref: 94261, follow these steps:

- 1-Separate the thigh and calf strips of foam of the orthosis, separating the articulations.
- 2-Position the thigh and calf strips, cutting off any excess and adjusting it on the leg, then following the same instructions as in point 3.

PRECAUTIONS

The constructive material is inflammable.

Do not expose these products to situations where they could ignite. If this were the case remove them quickly and use the necessary measures to put it out. The material used is hypoallergenic, however we cannot guarantee 100% that in some determined cases it will not produce allergies, if this occurs, remove the product and see the doctor that prescribed it.

In case of minor irritations due to sweat, we recommend the use of a cotton sock, to separate the skin from contact with the material. For small inconvenience caused by sweat, we recommend using an interface to separate the skin from contact with the tissue. If rash, irritation or swelling to remove the product and consult a doctor or prosthetist. Contraindicated in open scars with swelling, redness and accumulation heat.

RECOMMENDATION-WARNINGS

These products are to be used only by those which have been medically prescribed to do so. Persons unconnected to such prescriptions should not use them. When the product is not being used, store in its original packaging. In order to dispose of the packaging and the product, strictly comply with the legal norms of your community.

The use of the orthosis is conditioned to the recommendations of the prescribing doctor, therefore it must not be used for any other means than those prescribed.

For the validity of the guarantee, the issuing establishment should comply with this sheet of instructions.

MANUFACTURE-CHARACTERISTICS

The material features used in production have been tested and approved and meet all the European quality requirements.

All the products are produced with high quality materials and offer an unbeatable comfort and quality of use. All products offer restraint, stability and compression for the optimum treatment of the pathologies for which they have been designed.

RECOMMENDATIONS FOR PRESERVATION AND CLEANING

Stick the velcros to each other (if the orthosis has any), wash periodically by hand with a damp cloth and a neutral soap. To dry the orthosis, use a dry towel to absorb most of the dampness and leave to dry at room temperature. Do not hang up or iron and do not expose to direct heat sources such as stoves, heaters, radiators, direct sun light etc. During use or during washing do not use alcohols, ointments or dissolvent liquids. If the orthosis is not properly dried any detergent residues could irritate the skin and deteriorate the product.

GUARANTEES

ORLIMAN, S.L.U. Guarantees all its products, as long as these have not been manipulated or altered from their original state. Does not guarantee those products in which as a result of misuse, deficiencies or breakage of any kind, its original state has been altered. If you observe any deficiency or anomaly, communicate this immediately to the establishment from which it was obtained in order for it to be changed.

ORLIMAN, S.L.U. Appreciates your choice and wishes you a speedy recovery.

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