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## English

### INSTRUCTIONS FOR USE AND PRESERVATION

Dear customer,

Thank you for your trust in this Orliman product. You have acquired a quality product with a high medicinal value. Please, read the instructions carefully. If you have any doubt, get in contact with your doctor, your orthopaedist specialist or our customer service department. Orliman appreciates your choice and wishes you a speedy recovery.

#### REGULATION

These products satisfy the European Medical Device Directive 93/42/EEC (RD. 1591/2009). A Risk Analysis has been carried out (UNE EN ISO 14971) minimising all the existing risks. The trials have been carried out in accordance with the European regulation UNE-EN ISO 22523 of Prosthesis and Orthosis.

#### SUITABLE FOR

Protection and control following traumatism or surgical intervention, sprains and ligament injuries, as a method to replace plaster in the treatment of preulcerous or ulcerous conditions of the foot's plantar surface. Acute plantar fasciitis and in those cases which require a uniform distribution of body weight on the plantar surface and immobilisation of the tibia-tarsus articulation.

#### FITTING INSTRUCTIONS

**THE PRODUCT SHOULD BE ADJUSTED BY YOUR ORTHOPAEDIC TECHNICIAN OR IN ACCORDANCE WITH DOCTOR'S ORDERS.**

In order to obtain a greater degree of therapeutic effectiveness in the different pathologies and to prolong the life and utility of the product, it is fundamental that the correct size adequate for each patient or user is chosen (enclosed in the packet is a diagram with all the sizes and the equivalent in centimetres). An excess compression can produce intolerance, therefore we recommend that the compression be regulated to an optimum degree.

#### A - Follow the fitting instructions below:

It is advisable to fit the walker with the patient seated and over a sock, preferably cotton.

- 1-Unfasten the Velcro to separate the rigid structure from the surrounding cover. Then fit the surrounding cover around the patient's foot and leg as tight as possible and fasten with the Velcro closure system. Fit the toe protection pad.
- 2-Next, insert the leg into the orthosis. Make sure that the foot makes contact with the base. (In the case of articulated walkers, make sure that the mediolateral splints are parallel to each other and centred over the load line of the leg and ankle). Fasten the rigid structure (or mediolateral splints, in the case of articulated walkers) to the surrounding cover with the Velcro straps located inside the structure.
- 3-Finally, close and tighten the orthosis with the Velcro straps (in the case of articulated walkers, fit the tibia and calf protection plates before tightening the straps).

#### B - Instructions for adjusting the joint (EST-083 and EST-086):

For proper functioning of the joint, it is advisable to not adjust the extension and flexion movements simultaneously.

- 1-Unlock the white regulator by sliding the orange piece out.
- 2-Then press the white piece at the same time as moving it to the desired angle. Listen for the click.
- 3-Slide the orange piece back to below the white piece to lock the position. This prevents the selected angle from moving or being altered.
- 4-Repeat the same process on the other side of the joint. To prevent damage to the mechanical joints and complications in the limb, the flexion and extension regulators should be set to the same positions.

(To lock the bars at a fixed angle of 90°, set the flexion and extension regulators to 0°).

#### C - Instructions for using the inflation system (EST-083 and EST-089):

- 1-To inflate, turn the valve in the INFLATE direction (clockwise). The valve must be tightly closed so that air does not leak.
- 2-Press the orange circular piece several times to reach the desired pressure.
- 3-To release air from the bag, turn the valve in the DEFLATE direction (anti-clockwise) until it feels comfortable.

### PRECAUTIONS

The constructive material is inflammable.

Do not expose these products to situations where they could ignite. If this were the case remove them quickly and use the necessary measures to put it out. The material used is hypoallergenic, however we cannot guarantee 100% that in some determined cases it will not produce allergies, if this occurs, remove the product and see the doctor that prescribed it.

In case of minor irritations due to sweat, we recommend the use of a cotton sock, to separate the skin from contact with the material. For small inconvenience caused by sweat, we recommend using an interface to separate the skin from contact with the tissue. If rash, irritation or swelling to remove the product and consult a doctor or prosthetist. Contraindicated in open scars with swelling, redness and accumulation heat.

### RECOMMENDATION-WARNINGS

These products are to be used only by those which have been medically prescribed to do so. Persons unconnected to such prescriptions should not use them. When the product is not being used, store in its original packaging. In order to dispose of the packaging and the product, strictly comply with the legal norms of your community.

The use of the orthosis is conditioned to the recommendations of the prescribing doctor, therefore it must not be used for any other means than those prescribed.

For the validity of the guarantee, the issuing establishment should comply with this sheet of instructions.


### MANUFACTURE-CHARACTERISTICS

The material features used in production have been tested and approved and meet all the European quality requirements.

All the products are produced with high quality materials and offer an unbeatable comfort and quality of use. All products offer restraint, stability and compression for the optimum treatment of the pathologies for which they have been designed.

### RECOMMENDATIONS FOR PRESERVATION AND CLEANING

Stick the velcros to each other (if the orthosis has any), wash periodically by hand with a damp cloth and a neutral soap. To dry the orthosis, use a dry towel to absorb most of the dampness and leave to dry at room temperature. Do not hang up or iron and do not expose to direct heat sources such as stoves, heaters, radiators, direct sun light etc. During use or during washing do not use alcohols, ointments or dissolvent liquids. If the orthosis is not properly dried any detergent residues could irritate the skin and deteriorate the product.

### GUARANTEES

ORLIMAN, S.L.U. Guarantees all its products, as long as these have not been manipulated or altered from their original state. Does not guarantee those products in which as a result of misuse, deficiencies or breakage of any kind, its original state has been altered. If you observe any deficiency or anomaly, communicate this immediately to the establishment from which it was obtained in order for it to be changed.

ORLIMAN, S.L.U. Appreciates your choice and wishes you a speedy recovery.

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