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English**INSTRUCTIONS FOR USE AND PRESERVATION**

Dear customer,

Thank you for your trust in this Orliman product. You have acquired a quality product with a high medicinal value. Please, read the instructions carefully. If you have any doubt, get in contact with your doctor, your orthopaedist specialist or our customer service department. Orliman appreciates your choice and wishes you a speedy recovery.

REGULATION

These products satisfy the European Medical Device Directive 93/42/CEE (RD. 1591/2009). A Risk Analysis has been carried out (UNE EN ISO 14971) minimising all the existing risks. The trials have been carried out in accordance with the European regulation UNE-EN ISO 22523 of Prosthesis and Orthosis.

SUITABLE FOR

Specially indicated in patients affected by kyphotic attitudes, vertebral wedging, dorsolumbar traumatism, infectious illnesses which affect the dorsal rachides, slipped disks, postoperative, osteoporosis and in all those cases in which the object is the correction of the kyphotic curve, and the antepulsion of shoulders.

Some lower dorsalgias T-8, T-12, associated with lumbalgias, require temporary immobilization (movement of the patient, domestic and work-related mechanical efforts, etc.)

In those cases in which they partially need thoracic contention at some time of the day, depending on the activity.

The dorsal module can always adapt upon an Evotec lumbrosacral support.

FITTING INSTRUCTIONS

THE PRODUCT SHOULD BE ADJUSTED BY YOUR ORTHOPAEDIC TECHNICIAN OR IN ACCORDANCE WITH DOCTOR'S ORDERS.

In order to obtain a greater degree of therapeutic effectiveness in the different pathologies and to prolong the life and utility of the product, it is fundamental that the correct size adequate for each patient or user is chosen (enclosed in the packet is a diagram with all the sizes and the equivalent in centimetres). An excess compression can produce intolerance, therefore we recommend that the compression be regulated to an optimum degree.

To fit observe the following aspects:

Putting on the Dorso Evotec :

The dorsal piece is designed to be put on in combination with the Evotec high and low girdle.

Once the model and the correct size for the patient are selected, we can assemble the dorsal piece with the Evotec girdle. Prior to this, the lumbar pad must be removed from the girdle, as it is not possible to put it on together with the Dorso Evotec.

1-2: Place the Evotec girdle upon a smooth surface and open it completely; then stick on the dorsal piece, with the internal velour of the evotec and the microhook of the Dorso Evotec, making sure that the 2 central stays of the girdle coincide with the 2 central stays of the dorsal piece, always starting to join from the lower part of the lumbar girdle; the height can be regulated from a minimum of 47 cm to a maximum of 53 cm in the low model and from a minimum of 47 cm to a maximum of 57 cm in the high model. (see drawing).

This will allow us to have a contact and adherence area that is large enough for both pieces to be stable together. If we wish to increase the height, the area of adherence will be minimal, losing solidity in the combination of both components.

3-4-5: Once both pieces are adhered at the desired height, adapt it to the patient, starting by fastening the abdominal microhook lock until the desired compression is obtained.

Then, pull the side lumboabdominal tensors until obtaining the desired compression, fastening them to the front of the girdle with the microhook lock.

6: Pull the straps over the shoulders and under the axilla, crossing them over the back and up to the abdomen, joining them together with the lock system.

7-8: Regulation of the length of the straps:

At its end, together with the lock ring, it has an adhesive system which, when removed, will leave the end of the cinch free; the excess length can then be cut (see drawing).

Once this operation is carried out, join the end again to the adhesive system.

9: Last, adjust the lock, joining both rotational buckles, and stretching until hearing a click. The adjustment of this lock must be on the central front part of the Evotec girdle.

Adaptation of the module in Evotec Forte thermoplastic with the Dorso Evotec dorsal piece:

A: We will first joint the 2 side-fastening cinches to the Dorsal piece.

B: We will then join the thermoplastic module to the dorsal piece, with the adhesive surfaces of both modules to fasten them laterally with the 2 microhook fastening cinches attached to the thermoplastic module.

The height may be regulated from a minimum of 47 cm to a maximum of 53 cm in the low model and of 47 cm to a maximum of 57 cm in the high model.

C: Then adapt the Evotec girdle, clasping both components (thermoplastic module and dorsal module), thus proceeding to lock it, following the steps described previously in the placement of the Dorsal piece, from point no. 3.

PRECAUTIONS

The constructive material is inflammable.

Do not expose these products to situations where they could ignite. If this were the case remove them quickly and use the necessary measures to put it out. The material used is hypoallergenic, however we cannot guarantee 100% that in some determined cases it will not produce allergies, if this occurs, remove the product and see the doctor that prescribed it.

In case of minor irritations due to sweat, we recommend the use of a cotton sock, to separate the skin from contact with the material. For small inconvenience caused by sweat, we recommend using an interface to separate the skin from contact with the tissue. If rash, irritation or swelling to remove the product and consult a doctor or prosthetist. Contraindicated in open scars with swelling, redness and accumulation heat.

RECOMMENDATION-WARNINGS

These products are to be used only by those which have been medically prescribed to do so. Persons unconnected to such prescriptions should not use them. When the product is not being used, store in its original packaging. In order to dispose of the packaging and the product, strictly comply with the legal norms of your community.

The use of the orthosis is conditioned to the recommendations of the prescribing doctor, therefore it must not be used for any other means than those prescribed.

For the validity of the guarantee, the issuing establishment should comply with this sheet of instructions.

MANUFACTURE-CHARACTERISTICS

The material features used in production have been tested and approved and meet all the European quality requirements.

All the products are produced with high quality materials and offer an unbeatable comfort and quality of use. All products offer restraint, stability and compression for the optimum treatment of the pathologies for which they have been designed.


RECOMMENDATIONS FOR PRESERVATION AND CLEANING

Stick the velcros to each other (if the orthosis has any), wash periodically by hand with a damp cloth and a neutral soap. To dry the orthosis, use a dry towel to absorb most of the dampness and leave to dry at room temperature. Do not hang up or iron and do not expose to direct heat sources such as stoves, heaters, radiators, direct sun light etc. During use or during washing do not use alcohols, ointments or dissolvent liquids. If the orthosis is not properly dried any detergent residues could irritate the skin and deteriorate the product.

**GUARANTEES**

ORLIMAN, S.L.U. Guarantees all its products, as long as these have not been manipulated or altered from their original state. Does not guarantee those products in which as a result of misuse, deficiencies or breakage of any kind, its original state has been altered. If you observe any deficiency or anomaly, communicate this immediately to the establishment from which it was obtained in order for it to be changed.

ORLIMAN, S.L.U. Appreciates your choice and wishes you a speedy recovery.

- For proper use of the orthosis it is advisable not to modify its height position once the Velcro straps are adjusted, since the fabric could tear. If you need to adjust the height of the orthosis, first loosen the Velcro straps, re-position the orthosis at the desired height, and adjust the tension of the Velcro straps.

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