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**English****INSTRUCTIONS FOR USE AND PRESERVATION**

Dear customer,

Thank you for your trust in this Orliman product. You have acquired a quality product with a high medicinal value. Please, read the instructions carefully. If you have any doubt, get in contact with your doctor, your orthopaedist specialist or our customer service department. Orliman appreciates your choice and wishes you a speedy recovery.

**REGULATION**

These products satisfy the European Medical Device Directive 93/42/CEE (RD. 1591/2009). A Risk Analysis has been carried out (UNE EN ISO 14971) minimising all the existing risks. The trials have been carried out in accordance with the European regulation UNE-EN ISO 22523 of Prosthesis and Orthesis.

**SUITABLE FOR**

Disc hernias. Osteoporosis. Spondyloarthritis. Spondyloarthritis with or without stability. Traumatism. Compression fractures. Degenerative disc disease. Post-operation fusion. Chronic muscular disease. Vertebral metastasis treatment. TlSO model also suitable for the correction of hyperkyphosis and the prevention of channel stenosis in bedridden patients. The Lsoc and Tlso models are particularly suitable for patients who are most of the time in sedestation.

**FITTING INSTRUCTIONS**

**THE PRODUCT SHOULD BE ADJUSTED BY YOUR ORTHOPAEDIC TECHNICIAN OR IN ACCORDANCE WITH DOCTOR'S ORDERS.**

In order to obtain a greater degree of therapeutic effectiveness in the different pathologies and to prolong the life and utility of the product, it is fundamental that the correct size adequate for each patient or user is chosen (enclosed in the packet is a diagram with all the sizes and the equivalent in centimetres). An excess compression can produce intolerance, therefore we recommend that the compression be regulated to an optimum degree.

**To fit observe the following aspects:**

With the patient in supine position.

**A-Sacrolumbar corset positioning:**

1. Introduce the Velcro strap through the fasteners without causing traction.
- 2-3. Position the orthosis on the patient and seal using the zip located on the right side.
- 4-5. First tighten the lower and then the upper regulation straps, until the required compression is obtained.
6. In this manner, the orthosis is perfectly adjusted to the patient.

**B-Thoracic corset adaptation.**

Adaptation of the straps for the retropulsion of the shoulders.

The thoracic-sacrolumbar orthosis is fitted with shoulder straps, the function of which is to provide shoulder retropulsion. The straps start at the intercostal area between the scapula and the armpit. When positioning, pass it over the back and under the armpit, thus supporting the shoulder on the front and upper section, and then return to the back, introducing it through the fastener located in the rear upper section of the orthosis. With the orthosis suspended from the shoulders, adapt the lower part of the lumbar thoracic orthosis as follows (see sacrolumbar corset positioning texts).

**Adjusting the straps:**

Once the lumbar area of the orthosis is adapted, adjust the shoulder straps which should criss-cross over the back (forming an X), at T-6 and T-8 level, displacing them towards the sides in search of the anterior abdominal wall. Then tighten the shoulder straps with both hands and use the micro-fastener to attach them to the Velcro which covers the abdominal wall. The straps have an end micro-fastener for adaptation to different morphologies; whenever the length is excessive, cut using the marks which have been designed for this purpose, thus removing any excess and preventing discomfort.

**C-Pendulum abdomen adaptation:**

The principal function of the pendulum abdomen corset is to elevate the abdominal wall. To do this, we recommend properly adapting the lower strap, as low as possible above the pubis. Fasten this and then tighten and fasten the upper strap. In some cases it may be necessary to adapt the orthosis with the patient in supine position with the knees in flexion, thus elevating the abdomen before positioning the orthosis.

**D-Adaptation plates:**

The orthosis has two thermoplastic plates, a lumbar plate and an abdominal plate, which should be introduced in the pockets located for this purpose on the rear face of the orthosis and on the abdominal skirt. The plates allow for small modifications or adaptations through moulding with a hot air pistol.

**Recommendations for the daily regulation and positioning of the orthosis:**

- 1-In order to remove the orthosis, first release the lower and upper attachment straps.
- 2-Then release the zip, completely opening the orthosis and allowing removal from the patient.
- 3-For subsequent positioning, position the orthosis on the patient and proceed to close the zip.
- 4-Proceed with the adjustment, first tightening the lower traction straps and then doing the same with the upper straps, until the required compression is obtained. The regulation should be made uniformly, pulling with both hands in an upward direction, ensuring there are no folds or creases.
- 5-Proceed to regulate as many times as necessary, ensuring that the straps maintain the same length.

**PRECAUTIONS**

The constructive material is inflammable.

Do not expose these products to situations where they could ignite. If this were the case remove them quickly and use the necessary measures to put it out. The material used is hypoallergenic, however we cannot guarantee 100% that in some determined cases it will not produce allergies, if this occurs, remove the product and see the doctor that prescribed it.

In case of minor irritations due to sweat, we recommend the use of a cotton sock, to separate the skin from contact with the material. For small inconvenience caused by sweat, we recommend using an interface to separate the skin from contact with the tissue. If rash, irritation or swelling to remove the product and consult a doctor or prosthetist. Contraindicated in open scars with swelling, redness and accumulation heat.

**RECOMMENDATION-WARNINGS**

These products are to be used only by those which have been medically prescribed to do so. Persons unconnected to such prescriptions should not use them. When the product is not being used, store in its original packaging. In order to dispose of the packaging and the product, strictly comply with the legal norms of your community.

The use of the orthosis is conditioned to the recommendations of the prescribing doctor, therefore it must not be used for any other means than those prescribed.

For the validity of the guarantee, the issuing establishment should comply with this sheet of instructions.

**MANUFACTURE-CHARACTERISTICS**

The material features used in production have been tested and approved and meet all the European quality requirements.

All the products are produced with high quality materials and offer an unbeatable comfort and quality of use. All products offer restraint, stability and compression for the optimum treatment of the pathologies for which they have been designed.

**RECOMMENDATIONS FOR PRESERVATION AND CLEANING**

Stick the velcros to each other (if the orthosis has any), wash periodically by hand with a damp cloth and a neutral soap. To dry the orthosis, use a dry towel to absorb most of the dampness and leave to dry at room temperature. Do not hang up or iron and do not expose to direct heat sources such as stoves, heaters, radiators, direct sun light etc. During use or during washing do not use alcohols, ointments or dissolvent liquids. If the orthosis is not properly dried any detergent residues could irritate the skin and deteriorate the product.

**GUARANTEES**

ORLIMAN, S.L.U. Guarantees all its products, as long as these have not been manipulated or altered from their original state. Does not guarantee those products in which as a result of misuse, deficiencies or breakage of any kind, its original state has been altered. If you observe any deficiency or anomaly, communicate this immediately to the establishment from which it was obtained in order for it to be changed.

ORLIMAN, S.L.U. Appreciates your choice and wishes you a speedy recovery.

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