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**English****INSTRUCTIONS FOR USE AND PRESERVATION**

Dear customer,

Thank you for your trust in this Orliman product. You have acquired a quality product with a high medicinal value. Please, read the instructions carefully. If you have any doubt, get in contact with your doctor, your orthopaedist specialist or our customer service department. Orliman appreciates your choice and wishes you a speedy recovery.

**REGULATION**

These products satisfy the European Medical Device Directive 93/42/CEE (RD. 1591/2009). A Risk Analysis has been carried out (UNE EN ISO 14971) minimising all the existing risks. The trials have been carried out in accordance with the European regulation UNE-EN ISO 22523 of Prosthesis and Orthosis.

**SUITABLE FOR****TN-260/TN-261**

Indications: Tendonitis, arthrosis and arthritis, support and protection for sports and work activities, slight contusions, minor trauma and instability. Inflammatory processes.

**TN-262**

Indications: Prevention of sporting and work injuries, trapeziometacarpal joint irritation associated with wrist injuries.

**TN-230**

Indications: Lateral and medial epicondylitis, tendonitis, support and protection for sports and work activities, slight contusions, minor trauma and instability. Inflammatory processes.

**TN-210**

Indications: Arthritis and arthrosis, support and protection for sports and work activities, slight contusions, minor trauma and instability. Inflammatory processes.

**TN-211**

Indications: Arthritis and arthrosis brace, support and protection for sports and work activities, slight contusions, minor trauma and instability. Inflammatory processes.

**TN-240/TN-241**

Indications: Slight contusions and sprains, tendonitis, support and protection for sports and work activities, minor trauma and instability. Inflammatory processes.

**FITTING INSTRUCTIONS**

**THE PRODUCT SHOULD BE ADJUSTED BY YOUR ORTHOPAEDIC TECHNICIAN OR IN ACCORDANCE WITH DOCTOR'S ORDERS.**

In order to obtain a greater degree of therapeutic effectiveness in the different pathologies and to prolong the life and utility of the product, it is fundamental that the correct size adequate for each patient or user is chosen (enclosed in the packet is a diagram with all the sizes and the equivalent in centimetres). An excess compression can produce intolerance, therefore we recommend that the compression be regulated to an optimum degree.

**To fit observe the following aspects:****TN-260**

Place the wrist support on your hand and introduce your thumb through the opening. Adjust it, thereby avoiding creases and folds. Excessive compression can cause interruptions to the blood flow.

**TN-261**

Open the first and second part of the wrist strap.

Insert your hand through the elastic ring. Adjust the padding to the palm of the wrist. Tie around the wrist, and fasten with the small hooks. Apply the desired tension by completely surrounding the wrist with the external elastic strap, and fasten with the small hook.

**TN-262**

Insert the thumb into the thumb guard, ensuring the seam coincides with the interdigital commissure.

Wrap the elastic wrist band around the wrist to obtain the desired tension and compression and fasten with the hook.

**TN-230**

In order to put it on correctly, introduce your arm through the widest part of the elbow support. Position it at the elbow joint, with the pad at the back of the joint. Excessive compression can cause interruptions to the blood flow.

**TN-210**

In order to put it on correctly, slide on the knee support, starting with the widest end, and place it over the knee joint by positioning the pad over the patella. Excessive compression can cause interruptions to the blood flow.

**TN-211**

In order to put it on correctly, slide on the knee support, starting with the widest end, and place it over the knee joint by positioning the pad over the patella, thereby ensuring that the braces are parallel to each other and that they coincide with the medial line of the leg. Excessive compression can cause interruptions to the blood flow.

To fit observe the following aspects, place the product in the corresponding joint segment, elbow, ankle, hand, knee and arm so that they are aligned or fitted to the muscular movement sequences. Adjust the products so that they are adapted to the natural configuration and don't produce irritation or discomfort during their use.

**TN-240**

In order to put the support on correctly, introduce your foot in the ankle support and position your heel in the area designed for this purpose so that the pad is hugging your Achilles tendon and the malleoli. Tighten it and adjust it, thereby avoiding creases and folds. Excessive compression can cause interruptions to the blood flow. Be sure that the ankle support is in the right direction. Use the seam at the ends as a reference, with the widest opening at the top.

**TN-241\* Stabilizing bandage****TN-241\*\* Antivaro stabilizing bandage tape**

The brace TN-241, can be adapted interchangeably for both left and right foot to a foot switch left and right, to change from one foot to the other, flip the bandage.

Slide your foot into the ankle just behind the metatarsal phalangeal joint.

Once positioned correctly the elastic band of the foot, adjust the elastic strap around the ankle embraced as shown drawings 2 and 3.

Once the desired tension of the elastic strap attach the velcro on the bandage.

**PRECAUTIONS**

The constructive material is inflammable.

Do not expose these products to situations where they could ignite. If this were the case remove them quickly and use the necessary measures to put it out. The material used is hypoallergenic, however we cannot guarantee 100% that in some determined cases it will not produce allergies, if this occurs, remove the product and see the doctor that prescribed it.

In case of minor irritations due to sweat, we recommend the use of a cotton sock, to separate the skin from contact with the material. For small inconvenience caused by sweat, we recommend using an interface to separate the skin from contact with the tissue. If rash, irritation or swelling to remove the product and consult a doctor or prosthetist. Contraindicated in open scars with swelling, redness and accumulation heat.

**RECOMMENDATION- WARNINGS**

These products are to be used only by those which have been medically prescribed to do so. Persons unconnected to such prescriptions should not use them. When the product is not being used, store in its original packaging. In order to dispose of the packaging and the product, strictly comply with the legal norms of your community.

The use of the orthosis is conditioned to the recommendations of the prescribing doctor, therefore it must not be used for any other means than those prescribed.

For the validity of the guarantee, the issuing establishment should comply with this sheet of instructions.

**MANUFACTURE - CHARACTERISTICS**


The material features used in production have been tested and approved and meet all the European quality requirements.

All the products are produced with high quality materials and offer an unbeatable comfort and quality of use. All products offer restraint, stability and compression for the optimum treatment of the pathologies for which they have been designed.

**RECOMMENDATIONS FOR PRESERVATION AND CLEANING**

Stick the velcros to each other (if the orthosis has any), wash periodically by hand with a damp cloth and a neutral soap.

To dry the orthosis, use a dry towel to absorb most of the dampness and leave to dry at room temperature. Do not hang up or iron and do not expose to direct heat sources such as stoves, heaters, radiators, direct sun light etc. During use or during washing do not use alcohols, ointments or dissolvent liquids. If the orthosis is not properly dried any detergent residues could irritate the skin and deteriorate the product.

**GUARANTEES**

ORLIMAN, S.L.U. Guarantees all its products, as long as these have not been manipulated or altered from their original state. Does not guarantee those products in which as a result of misuse, deficiencies or breakage of any kind, its original state has been altered. If you observe any deficiency or anomaly, communicate this immediately to the establishment from which it was obtained in order for it to be changed.

ORLIMAN, S.L.U. Appreciates your choice and wishes you a speedy recovery.

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